



Be Right™



Pocket Dissolved Oxygen (DO) Meter

หมายเลขผลิตภัณฑ์:

2968800

THB ราคา,
พร้อมจำหน่าย

ติดต่อ Hach

Pocket DO Meter measures dissolved oxygen to determine water quality in a variety of applications including wastewater, aquariums and fish hatcheries. Simultaneously displays DO and temperature readings. Unit holds and recalls minimum and maximum readings. Features automatic temperature compensation and calibration to air. Meter is ready for use and includes a soft carrying case with belt loop, two spare probe heads, electrolyte and 4 AAA batteries. Replacement electrolyte sold separately.

Sturdy design

Heavy duty polarographic probe measures dissolved oxygen and temperature.

ATC

Automatic temperature compensation readings.

Dual display

Displays D.O. and temperature at the same time.

Water resistant

IP-67 water resistant rating.

Max/Min

Maximum and minimum readings displayed.

ข้อมูลจำเพาะ

| | |
|-------------------------|--------------------------|
| ความละเอียด: | 1 mV |
| ความแม่นยำ 2: | ± 0.4 mg/L DO |
| ช่วง 2: | 0 - 20.0 mg/L DO |
| ช่วงอุณหภูมิในการทำงาน: | 0 - 50 °C |
| น้ำหนัก: | 0.176 kg |
| ปริมาณ: | 1 each |
| พลังงาน: | 4-AAA Alkaline Batteries |